

SPESIFIKASI TEKNIS DIRUI BF-6500 Automatic Hematology Analyzer

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	DIRUI BF-6500 Automatic Hematology Analyzer
2	Produsen	DIRUI INDUSTRIAL CO., LTD., P.R., CHINA
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 20205311164
5	Spesifikasi	
	Test item	WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, NEUT%, LYMPH%, MONO%, EO%, BASO%, NEUT#, LYMPH#, MONO#, EO#, BASO#, RDW-SD, RDW-CV, PDW, MPV, PCT, P-LCR
	Research parameter	BLAST#, IMM#, LEFT#, BLAST#, IMM%, LEFT%, ABNLYM#, NRBC#, ABNLYM%, NRBC%
	Test Principle	Semiconductor laser flow cytometry combined with cytochemical staining, impedance, enviromental friendly cyanide-free colorimetry
	Throughput	CBC+DIFF mode 60 samples/hour
	Analysis mode	CBC mode, CBC+DIFF mode
	Sample type	Whole blood, pre-dilution blood
	Sampling device	Manual sampling
	Data storage	With a storage capacity of 100.000 patient's result, including scattergram, histogram and patient's information. Multi QC methods including X, X-B, L-J, X and L-J totally 12 document, 400 results/document. X-B can save 1000 results.
	Language	English
	Display type	External computer
	Expansion function	USB port, internet port, support U-Disk, printer, mouse and keyboard, etc
	Working Condition	Temperature 18-30°C , humidity ≤ 75%
	Power	100-240VAC 50-60Hz