

TECHNICAL SPECIFICATION

Measuring range	0 to 30 mg/dl (total bilirubin)
Response	3 second
Correcting hemolysis	0 to 250 mg/dl
Concentration alarm	Indicates 999 for 31 mg/dl or more
Turbidity alarm	Turbidity of 250 mg/dl or more
Filters	455 nm & 575 nm
0-adjustment	Automatically adjusted by pushing a button
Sample blood volume	50 to 60 μ l (whole blood)
Cell	Designated capillary tube
Display	3-digit 7-segment red LED
Light source lamp	White LED
Dimensions	230 (W)x200(D)x100(H) mm
Weight	2.0 kg
Power supply	AC 90-240V 50/60Hz 15VA