

SPESIFIKASI TEKNIS CoaData 501

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	CoaData 501
2	Produsen	LABITEC (LABOR BIO MEDICAL TECHNOLOGY) GmbH, Germany
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 20205704283
5	Instrument type	Analyzer for determination of plasmic clotting
6	Measuring Principle	Turbodensitometric with automated start function
7	Programmed methods	PT, in sec. %, RATIO, INR aPTT, in sec, g/L, and mg/L, Thrombin Time in sec Intr. Factor, in % Extr. Factor, in %
8	Incubation block	Heated to 37°C (controlled at 37.4°C ± 0.3°C) 1 measuring cell with magnetic mixer motor 1 reagent position diameter 23.0 mm 4 cuvette incubation positions
9	Display	1 lines with 8 characters, liquid crystal display
10	Control panel	Membrane keypad
11	Light source	LED, light emitting diode
12	Test throughput	PT 30/h, aPTT 15/H, +/- 10 test/h
13	Calibration curve	Memory of 9 meas, points per assays
14	Dimensions	9,7 x 21,2 x 5,2 cm (WxDxH)
15	Memory	A maximum of 7 test methods can be store
16	Light protection cap	With yellow dot
17	Weight	0,6 kg
18	Power Requirement, VAC	9.6 VA ; 12 V DC
19	Interfaces	RS 232 C (optional printer) Chip CARD reader