

SPESIFIKASI TEKNIS Coatron A4

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	Coatron A4
2	Produsen	TECO GmbH, Germany
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 20205510828
5	Spesifikasi	
	Dimension	: 85 x 60 x 70 cm
	Weight	: 50 kg
	Power Input	: 90-240 Vac / 50-60 Hz
	Ambient condition	: 18-30°C, up to 80% r.H no direct sun emission CE Approved
	Safety	: EN 61010-1
	EMC	: EN 60601-1-2
	Calculate unit	: %, INR, R, mg/dL, µg/L
	Other Specials	: - Multilingual dialogue - Auto calibration - Autodilution (Fibrinogen, Factors) - Cuvettes for up to 160 determinations - 2 racks for up to 24 primary tubes - Emergency/STAT - QC controlled - Printout of instrument status - Printout of various service information - Bi-directional interface - Software updates via PC - Reaction curves
	Optional available	: - Trolley for the Coatron A4 - Benefits : enlargeds "Waste Management", Instrument is mobil, Keeps laboratory table free - Tecam; PC Software for data management