

SPESIFIKASI TEKNIS Coatron M4

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	Coatron M4
2	Produsen	TECO GmbH, Germany
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 20205805141
5	Spesifikasi	
	Dimensions	
	L x W x H	: 280 x 225 x 115mm
	Weight	: 2.30kg
6	Power Supply (External)	
	Input	: 100-240~, 47-63Hz
	Power	: 42W
	Approvals	: CE, EN60601-1, UL
7	Incubation Block	
	Temperature	: 37.0°C ± 0.2 °C
	Samples	: 12 Cuvettes
	Reagent	: 1 vial Ø30.0mm, stirred 1 vial Ø22.55mm 1 vial Ø16mm
8	Measurement	
	Optic	: 4 Channel with AutoSense
	Methods	: Turbidity, Kinetic, Endpoint
	Algorithm	: Based on optical density
	Processing	: 18Bit ADC + Digital Filter
	AutoSense	: LED intensity controlling
	Starting	: manual, optical, autopipette
9	Performance	
	Precision	CV <5%
	Throughput	PT ~150/h APTT ~50/h