

SPESIFIKASI TEKNIS DL-96 II

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	Zhuhai DL-96II Auto Microbial Identification and Antimicrobial Susceptibility Testing System and Accessories
2	Produsen	ZHUHAI DL BIOTECH CO., LTD., China
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 20301510019
5	Testing Principles	Identification: Colorimetry Susceptibility testing : Turbidimetry
6	Throughput	60-120 cards per hour
7	Optimal configuration	Automated sampling device, sampling 60-80 cards per hour
8	Printer	External printer
9	Interface	Internal PC, support LIS/HIS, WHONET
10	Working Environment	Environmental temperature : 50C-400C Relative humidity : 10%-90%
11	Atmospheric pressure	76kPa - 106kPa
12	Power supply	AC220V±22V, 50Hz±1Hz
13	Dimension & Weight	700mm x 493mm x 465mm ; 31.2 kg