

SPESIFIKASI TEKNIS Dimex Jr

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	Dimex Jr
2	Produsen	TECO GmbH, Germany
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 20207804744
5	Spesifikasi	
	Dimensions	
	L x W x H	: 245 x 130 x 65mm
	Weight	: 0.55kg
6	Power Supply (External)	
	Input	: 100-240~, 50/60Hz
	Ouput	: 12Vdc, 1,25A max
	Approvals	: CE, UL
7	Incubation Block	
	Temperature	: 37.0°C ± 0.2 °C
	Samples	: 6 Cuvettes
	Reagent	: 2 vial Ø11.0
8	Optic Block	
	One channel	~660nm
9	Measurement	
	Optic	: optical, with Auto-Sense
	Methods	: Turbidity, Kinetic, Endpoint
	Algorithm	: Based on optical density
	Processing	: 14Bit ADC + Digital Filter
	AutoSense	: LED intensity controlling
	Starting	: manual or optical
10	Performance	
	Precision	PT-B CV: < 10%
		DD CV: <12%
	Thoughtput	PT-B ~25/h
		DD ~25/h