

## **TECHNICAL SPECIFICATION**

<b>Photometer</b>	Photometric range : 0 to 3.0 OD Spectral range : 400nm to 700nm Read time : under 2 minutes per plate
<b>Incubator</b>	Temp range : up to 40 <sup>0</sup> C (104 <sup>0</sup> F)
<b>Pipetor</b>	Min/Max vol : 1 $\mu$ L to 300 $\mu$ L Precision (single dispense) : better than 10% Precision (multi dispense) : better than 10%
<b>Orbital Shaker</b>	Amplitude : 1mm Adjustable speed : 100-900 rpm Voltage : 24 volts Adjustable shaking time : 0 to 2 hrs
<b>Power</b>	100-265 V, 50-60 Hz
<b>Dimensions</b>	64 (W) x 57 (D) x 45 (H) cm
<b>Weight</b>	28 kg