

SPESIFIKASI TEKNIS K-Lite 3 Electrolyte Analyzer

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	K-Lite 3 Electrolyte Analyzer
2	Produsen	MEIZHOU CORNLEY HI-TECH Co., LTD, CHINA
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 20102011526
5	Features	Self made maintenance free electrodes Sleep mode to save reagents Rapid test speed ≤ 30 seconds/test Minimum volume sample consumption $\leq 65 \mu\text{l}$ Reagent Pack to avoid Pollution Online instruction provided for self problem solving ≥ 1000 record storage Automatic one point and two points calibration with additional manual calibration User friendly software
6	SAMPLE TYPE	Plasma, Serum, Whole Blood, Urine
7	TEST RANGE AND ACCURACY	K+ 0.50~10.00 mmol/L Na+ 20.00~200.00 mmol/L Cl- 20.00~200.00 mmol/L
8	Measuring Principle	Ion Selective Electrode (ISE)
9	Display	LCD Display
10	Computer interface	RS232
11	Test Result	1000 records
12	Printer	Thermal Printer
13	Dimension (L x W x H), cm (in)	30 x 26 x 36
14	WEIGHT, kg (lb)	7.6
15	POWER REQUIREMENTS, VAC	100 VAC Power supply 100-240 VAC at 50/60 Hz