

## SPESIFIKASI TEKNIS K-Lite 5 AU Electrolyte Analyzer

No.	Deskripsi	Uraian/Keterangan
1	<b>Nama Produk</b>	K-Lite 5 AU Electrolyte Analyzer and Accessories
2	<b>Produsen</b>	MEIZHOU CORNLEY HI-TECH Co., LTD, CHINA
3	<b>Distributor</b>	PT. Setia Anugrah Medika
4	<b>Nomor izin edar</b>	AKL 20101010042
5	<b>Features</b>	Self made maintenance free electrodes Sleep mode to save reagents Rapid test speed 80 test/hour Minimum volume sample consumption ≤ 65 µl Reagent Pack to avoid Pollution Online instruction provided for self problem solving ≥ 1000 record storage Automatic one point and two points calibration with additional manual calibration User friendly software Autosampler
6	<b>SAMPLE TYPE</b>	Plasma, Serum, Whole Blood, Urine
7	<b>TEST RANGE AND ACCURACY</b>	pH        6.00~9.00 K+        0.50~10.00 mmol/L Na+       20.00~200.00 mmol/L Cl-        20.00~200.00 mmol/L Ca++     0.30~5.00 mmol/L
8	<b>Measuring Principle</b>	Ion Selective Electrode (ISE)
9	<b>Display</b>	LCD Display
10	<b>Computer interface</b>	RS232
11	<b>Test Result</b>	1000 records
12	<b>Printer</b>	Thermal Printer
13	<b>Dimension (L x W x H), cm (in)</b>	30 x 26 x 36
14	<b>POWER REQUIREMENTS, VAC</b>	100 VAC Power supply 100-240 VAC at 50/60 Hz