

## SPESIFIKASI TEKNIS RAMP 200

No.	Deskripsi	Uraian/Keterangan
1	<b>Nama Produk</b>	RAMP 200 System
2	<b>Produsen</b>	RESPONSE BIOMEDICAL CORPORATION, Canada
3	<b>Distributor</b>	PT. Setia Anugrah Medika
4	<b>Nomor izin edar</b>	AKL 10102213337
5	<b>Spesifikasi</b>	
	<b>Test Module &amp; Control Module</b>	
	Berat	2.0 kg (4.4 lbs) for each Control Module and Test Module
	Dimensi	20 cm x 19.5 cm x 10.5 cm (pxlxt)
	Input Power Requirements	12V DC, 120W
	Power Supply	100-240 VAC, 50-60 Hz
	Keypad	Touch sensitive VGA display
	Display	640x480 pixel, 5,7-6,5 " diag (14.48~16.51 diag)
6	<b>Parameter Pemeriksaan</b>	Troponin I, CKMB (mesh), Myoglobin, Nt-proBNP, D-dimer, Influenza A/B, Respiratory Syncytial Virus (RSV), Anthrax, Botulinum Toxin, Ricin, Pox, West Nile Virus
7	<b>Prinsip Test</b>	Immunofluorescence (Fluorometer that analyzes results produced by immunoassay specified by Response Biomedical Corporation)
8	<b>Features</b>	High throughput , multi-port laboratory reader (18-36 test/hour with 3 TM)
		Modular, expandable and flexible design (expandable up to 3 TMs)
		self-diagnostic and internal quality control checks
		No calibration or maintenance required
		Run multiple test on the same or different patients simultaneously
		Stores up to 300 samples, 300 IQC and 300 LQC results