

SPESIFIKASI TEKNIS RAMP READER

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	RAMP Reader System
2	Produsen	RESPONSE BIOMEDICAL CORPORATION, Canada
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 10102213338
5	Spesifikasi	
	Berat	2.1Kg
	Dimensi	27cm x 25cm x 15 cm (pxlxt)
	Batere	Rechargeable - Nickel Cadmium (Ni-Cad) 12.5V, 1400mAh
	Power Requirement	100-240 VAC, 50-60 Hz Output power 15-18VDC, 1A max
	Layar	2 baris dengan 20 karakter huruf dan angka , Liquid crystal display (LCD)
	Keypad	10 tombol alphanumeric, 2 tombol fungsi spesial
	External Connection	RS-232, Mini DIN for Barcode wand
6	Parameter Pemeriksaan	Troponin I, CKMB (mesh), Myoglobin, Nt-proBNP, D-dimer, Anthrax, Botulinum Toxin, Ricin, Pox, West Nile Virus
7	Prinsip Test	Immunofluorescence (Fluorometer that analyzes results produced by immunoassay specified by Response Biomedical Corporation)
8	Features	Fully portable (battery and/or AC power) , lightweight reader, runs up to 100 test on one charge
		self-diagnostic and internal quality control checks
		No calibration or maintenance required
		Stores up to 500 samples, 65 IQC results and 35 LQC results
		Stores lot information - up to and including 50 lot
		Data can be exported to a computer or LIS