

SPESIFIKASI TEKNIS ThinPrep 2000 Processor

No.	Deskripsi	Uraian/Keterangan
1	Nama Produk	ThinPrep 2000 Processor
2	Produsen	HOLOGIC, INC., USA
3	Distributor	PT. Setia Anugrah Medika
4	Nomor izin edar	AKL 10203906427
5	Spesifikasi	
6	Intended use	ThinPrep 2000 Processor is an automated slide preparation unit used with the ThinPrep Test for Cervical cancer testing, and also be used for non-gynecological specimens such as for the diagnosis of cancers of the lung, bladder, and gastrointestinal tract, and the preparation of fine needle aspiration (FNA) of thyroid and breast
7	Throughput	25 samples per hour. Up to 50.000 samples per year
8	Dimension	46 cm x 38 cm x 50cm (pxlxt)
9	Weight	19 kg
10	Operating Temp	15-32°C
11	Electrical voltage	220/240 VAC at 1 amp
12	Power	Max 200 watts