

## SPESIFIKASI TEKNIS DL-Bt 64

No.	Deskripsi	Uraian/Keterangan
1	<b>Nama Produk</b>	ZHUHAI DL-BT 64 Auto Blood Culture Detection System and Accessories
2	<b>Produsen</b>	ZHUHAI DL BIOTECH CO., LTD., China
3	<b>Distributor</b>	PT. Setia Anugrah Medika
4	<b>Nomor izin edar</b>	AKL 10302413182
5	<b>Testing Principles</b>	Colorimetry & Carbonate sensor
6	<b>Sample Type</b>	Blood and various sterile body fluid
7	<b>Capacity</b>	64 bottles
8	<b>Detection Method</b>	<p>1. Adopt continuous swinging vibration mode increase bacterial growth, which enhances the detection rate of microorganism, even fastidious ones.</p> <p>2. Each vial cell has its own independent detector, which continuously monitors the real-time status of bottles.</p> <p>3. The instrument detects the bottles every 10 minutes to optimize detection time and increase its precision.</p>
9	<b>Bottle Loading</b>	Loading and unloading at any time
10	<b>Monitor</b>	Internal LCD 8 inch touch screen
11	<b>Working Environment</b>	<p>Environmental temperature :5<sup>0</sup>C - 40<sup>0</sup>c</p> <p>Relative humidity : 10% - 90%</p>
12	<b>Atmospheric pressure</b>	76 kPa - 106 kPa
13	<b>Power supply</b>	AC220V ± 22V, 50Hz ± 1 Hz
14	<b>Dimension &amp; Weight</b>	575 x 560 x 950 (mm); 80 kg